

TRANSCRANIAL ECOCOLOR DOPPLER from methodology to clinical applications


Transcranial Ecocolor Doppler is currently the elective Imaging method for noninvasive study of intracranial arteries and their hemodynamics. Feasible in adults and children, it is indicated for the early diagnosis of acute brain disease and in the follow-up of chronic cerebrovascular states.

TCD: EXAMINATION TECHNIQUE

With the patient in the supine position and with the operator behind the patient's head, it is possible to visualize the intracranial vessels through the 3 acoustic windows, areas of the skull, made up of the bones that allow the penetration of the US:

Ultrasound windows: Temporal, Occipital, Frontotemporal

Among them, the most interesting for diagnostic use is the temporal window (Fig. 1,2).





Rectangular area on the squamous part of the temporal bone subdivided in 3 zones: front, middle and rear.

Middle cerebral artery (MCA)

It collects about 60-70% of the blood coming from the internal carotid arteries (ICA), allowing, therefore, the estimate of blood flow to the hemisphere ipsilateral (Fig. 3).


It allows to explore (Fig. 4):



Middle cerebral arteries (MCA); in the M1 and M2 Anterior cerebral arteries (ACA) sections; in A1 and A2, Posterior cerebral arteries (PCA), P1 and P2.

In most patients, identification is possible through the 'butterfly wings sign' (Fig. 5).

In intracranial vessels, the Dopplersignal highlights a flow with a predominant diastolic component (Fig. 6).



TCD: INDICATIONS

From the spectral analysis of the Doppler signal it is possible to calculate the cerebral flow rate. Cerebral (Mean CBPF) (Table 1, 2)

Table 1 - Factors influencing cerebral blood flow velocity

Factor	Change in CBPF
Age	Increase up 6-10 yr than decrease
Sex	Woman > man
Pregnancy	Decrement in the III trimester
Hematocrit	Increase with decreasing Hct
PCO ₂	Increase with increasing PCO ₂
MAP	Atrial pressure increase with increasing MAP

CBF: Central blood flow velocity; MAP: Mean arterial pressure.

Table 2 - Mean cerebral blood flow (cm/s) related to age

Artery	Age 20-40 yr	Age 40-60 yr	Age > 60 yr
Anterior cerebral artery	56-80	53-61	44-51
Middle cerebral artery	74-91	72-73	58-59
Posterior cerebral artery P1	43-57	41-56	37-47
Posterior cerebral artery P2	43-51	40-57	37-47
Vertebral artery	37-51	29-50	30-37
Basilar artery	39-58	27-56	29-47

$$\text{Mean CBPF} = (\text{PSV} + [\text{EDV} \times 2]) / 3$$

PSV = peak systolic velocity

EDV = end-diastolic blood flow velocity

Change of CBPF mean, may suggest vasospasm

Lindegaard Ratio (LR): (Table 3)

Table 3 - Intracranial arteries: severity of vasospasm

	WFV	LR unmodified
MCA or ICA vaso spasm (%)		
Mild (< 25)	120-140	3-6
Moderate (25-50)	15-190	3-6
Severe (> 50)	> 200	> 6
BA vaso spasm (%)		
Possible vaso spasm	70-85	2-2.49
Moderate (25-50)	> 35	2.5-2.99
Severe (> 50)	> 85	> 3

MCA: Middle cerebral artery; ICA: Internal carotid artery; LR: Lindegaard ratio;

BA: Basilar artery; WFV: Mean flow velocity; MCA: Middle cerebral artery; ICA: Internal carotid artery; LR: Lindegaard ratio; BA: Basilar artery; MPA: Mean flow index.

MCA mean CBPF/extracranial ICA mean CBPF (v.n. <3)

BA mean CBPF/left or right extracranial VA mean CBPF (v.n. <2)



TCD: CRYPTOGENIC STROKE AND t-RLS

In subjects aged < 55 years, in 40% of cases stroke is of cryptogenic origin. In these pts, the prevalence of Patent Foramen Ovale (PFO), responsible for paradoxical embolism due to transient Right to Left shunt (t-RLS), is about 40%.

Examination technique (Fig. 7)


t-RLS is highlighted through shaken saline solution that behaves as a contrast.

Contrast agent: 9 mL of normal saline solution with 1 mL of air or blood, shaken up about 10 times through a system constituted by two 10 mL syringes linked by a 3-way stopcock. The agitated solution is administrated into the antecubital vein by an 18-gauge. The patient is then invited to perform a forced expiration against the closed glottis for a minimum of 10 seconds.

Table 4 - Grade of transient right to left shunting based on microembolic signals grading score

Grade transient shunt	MES
No shunt	0
Low grade shunt	1-10
Moderate grade shunt	11-25
High grade shunt	> 25 (showed) or uncountable (curtain effect)

C-TCD, placing the PW Doppler slider on the MCA, show the passage of 'Microembolic Signals-MES', indirect manifestations of t-RLS.



Content induced transcranial Doppler as a first-line screening test in the setting of a suspected transient ischemic attack (TIA). Transcranial Doppler (TCD) is a non-invasive technique for measuring blood flow in the brain. It can be used to detect changes in blood flow in the major arteries and veins of the brain.

Transcranial Doppler Ultrasound: A non-invasive and safe method for assessing blood flow in the brain. It uses high-frequency sound waves to measure the speed and direction of blood flow in the major arteries and veins of the brain. Transcranial Doppler (TCD) is a non-invasive technique for measuring blood flow in the brain. It can be used to detect changes in blood flow in the major arteries and veins of the brain.

esaote