



Rev. 19

March 2020

## **PROBES AND CONSUMABLES**

*141003500*

# Appendix C - Probes and Agents

## Agents and Systems

The following tables list the **MyLab** probes and the recommended cleaning, disinfection and sterilization agents and systems.

The tables show both chemical (and physical) compatibility and biological effectiveness of the method.

**NOTE** *The listed products may not be available in all regions. For probe availability, please contact your local representative; for cleaning, disinfectant or gel availability, contact the cleaning agent or gel manufacturer.*

OK means the use of the agent/system is permitted.

N means the use of the agent/system is not permitted.

- means the agent/system has not been tested.

**CAUTION**

**Do not use not tested agents/systems either to clean or to disinfect or to sterilize the probe. The probe could be damaged if such agents are used.**

Cleaning, disinfection and sterilization are allowed up to the immersible part limit.

**NOTE** *Any damage caused by the use of not recommended agents/systems is not covered by the warranty.*

**WARNING**

**The listed disinfection/sterilization agents and systems are recommended because of chemical compatibility with the probe materials.**

The tables in the next pages indicate which agents/systems have been positively tested for chemical compatibility only. Follow the guidelines and

**recommendations of the manufacturer for the biological effectiveness of agents/systems.**

**Use of solutions other than those referenced is not recommended. They may damage the probe housing or the acoustic lens. Esaote takes no responsibility for damage caused by using non-approved products.**

**Follow the instructions provided by the manufacturer of the agent for proper use. Observe specifically soak times and dilution rates.**

**Overexposure to the disinfection fluid can damage the probe.**

**Personnel should adopt all necessary protective measures during the probe cleaning, disinfection and sterilization processes (for example gloves, protective glasses).**

**Never attempt to clean or disinfect the probes while they are connected to the system.**

The probes have resulted compatible, in terms of chemical aggression, with reprocessing agents/methods listed in the tables below for average use of at least 7 years corresponding to not less than 4000 cycles. For sterilization, compatibility test evaluation is limited to 100 cycles.

In any case reprocessing methods are harsh and can shorten the life of the equipment owing to many unpredictable factors.

Equipment that undergoes reprocessing should be checked regularly.

Table C-1: Cleaning Agents

	Accel RTU wipes	Asept-wipe II	Cidezyme	CidezymeXTRA	Cleanisept-wipes forte	Enzol	Hexantos C+R	Klenzyme	Metrizyme	Microzid AF wipes	Microzid PAA wipes	Microzid Sensitive wipes	Mild soap	Percept RTU cleaner	Prolystica	Protex wipes	Sani-Cloth Active wipes	SaniZide plus spray	SaniZide plus wipes	Trionic D wipes	Trionic D spray	Tristel sporicidal wipes	Virusolve + wipes
2 MHz CW	OK	-	-	-	-	-	-	-	-	-	OK	-	OK	OK	-	-	-	-	-	-	-	-	-
5 MHz CW	OK	-	-	-	-	-	-	-	-	-	OK	-	OK	OK	-	-	-	-	-	-	-	-	-
AC2541	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
AL2442	OK	OK	-	-	OK	-	OK	-	-	-	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	-	OK
BC441	OK	OK	OK	OK	OK	OK	OK	-	OK	-	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	-	OK
BE1123	OK	-	-	-	-	-	-	-	-	-	OK	-	OK	OK	-	-	-	-	-	-	-	-	-
BL433	-	-	-	-	-	-	-	-	-	-	-	-	OK	-	-	-	-	-	-	-	-	-	-
C 1-8	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
C 2-9	OK	OK	-	-	OK	-	OK	-	-	-	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	-	OK
CA123	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
CA430E	OK	-	OK	OK	OK	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
CA431	OK	OK	OK	OK	OK	OK	OK		OK	-	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	-	OK
CA541	OK	-	-	-	-	-	-	-	-	-	OK	-	OK	OK	-	-	-	-	-	-	-	-	-
CA631	OK	OK	OK	OK	OK	OK	OK	-	OK	-	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	-	OK
E 3-12 <sup>a</sup>	OK	OK	-	-	OK	-	OK	-	-	-	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	-	OK
EC1123 <sup>a</sup>	OK	OK	-	-	OK	-	OK	-	-	-	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	OK	OK
EC123 <sup>a</sup>	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
HF CW	OK	-	-	-	-	-	-	-	-	-	OK	-	OK	OK	-	-	-	-	-	-	-	-	-
IH 6-18	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
IL 4-13	-	N	OK	OK	-	OK	-	OK	OK	-	-	-	OK	-	OK	-	-	-	-	-	-	OK	-
IOE323	-	N	OK	OK	-	OK	-	OK	OK	-	-	-	OK	-	OK	-	-	-	-	-	-	OK	-
IOT342	OK	OK	-	-	OK	-	OK	-	-	-	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	-	OK

	Accel RTU wipes	Asepti-wipe II	Cidezyme	CidezymeXTRA	Cleanisept-wipes forte	Enzol	Hexanios G+R	Klenzyme	Metrizyme	Microzid AF wipes	Microzid PAA wipes	Microzid Sensitive wipes	Mild soap	Percept RTU cleaner	Prolystica	Protex wipes	Sani-Cloth Active wipes	Santizide plus spray	Santizide plus wipes	Trionic D wipes	Trionic D spray	Tristel sporicidal wipes	Virusolve + wipes
L 3-11	OK	OK	-	-	OK	-	OK	-	-	-	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	-	OK
L 4-15	OK	OK	-	-	OK	-	OK	-	-	-	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	-	OK
L 8-24	-	-	-	-	-	-	-	-	-	-	-	-	OK	-	-	-	-	-	-	-	-	-	-
LA332E	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
LA435	-	OK	-	-	OK	-	OK	-	-	-	-	OK	OK	-	-	OK	OK	OK	OK	OK	OK	-	OK
LA523	-	OK	N	N	OK	N	OK	N	N	-	-	OK	OK	-	N	OK	OK	OK	OK	OK	OK	-	OK
LA533	OK	OK	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK		OK
LP323	-	N	OK	OK	-	OK	-	OK	OK	-	-	-	OK	-	OK	-	-	-	-	-	-	OK	-
LP 4-13	-	N	OK	OK	-	OK	-	OK	OK	-	-	-	OK	-	OK	-	-	-	-	-	-	OK	-
mC 3-11	OK	OK	-	-	OK	-	OK	-	-	-	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	OK	OK
P 1-5	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
P 2-9	OK	OK	-	-	OK	-	OK	-	-	-	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	-	OK
P 3-11	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
P2 3-11	OK	OK	-	-	OK	-	OK	-	-	-	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	-	OK
P2 5-13	OK	-	OK	OK	OK	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
PA023E	OK	-	OK	OK	OK	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
PA121E	OK	-	-	OK	OK	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
PA122E	OK	-	OK	OK	OK	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
PA230E	OK	-	OK	OK	OK	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
PA240	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
PA250	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
S2MCW	OK	-	-	-	-	-	-	-	-	-	OK	-	OK	OK	-	-	-	-	-	-	-	-	-
S2MPW	OK	-	-	-	-	-	-	-	-	-	OK	-	OK	OK	-	-	-	-	-	-	-	-	-
S5MCW	OK	-	-	-	-	-	-	-	-	-	OK	-	OK	OK	-	-	-	-	-	-	-	-	-

	Accel RTU wipes	Asepti-wipe II	Cidezyme	CidezymeXTRA	Cleanisept-wipes forte	Enzol	Hexanios G+R	Klenzyme	Metrizyme	Microzid AF wipes	Microzid PAA wipes	Microzid Sensitive wipes	Mild soap	Percept RTU cleaner	Polystica	Protex wipes	Sani-Cloth Active wipes	Santizide plus spray	Santizide plus wipes	Trionic D wipes	Trionic D spray	Tristel sporicidal wipes	Virusolve + wipes
SB2C41	OK	OK	OK	OK	OK	OK	OK	-	OK	-	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	-	OK
SB3123	OK	-	-	-	-	-	-	-	-	-	OK	-	OK	OK	-	-	-	-	-	-	-	-	-
SC3123	OK	OK	-	-	OK	-	OK	-	-	-	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	OK	OK
SC3421	OK	OK	OK	OK	OK	OK	OK	-	OK	-	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	OK	OK
SE3123 <sup>a</sup>	OK	OK	-	-	OK	-	OK	-	-	-	N	OK	OK	N	-	OK	OK	OK	OK	OK	OK	OK	OK
SE3133 <sup>a</sup>	OK	OK	-	-	OK	-	OK	-	-	-	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	OK	OK
SHFCW	OK	-	-	-	-	-	-	-	-	-	OK	-	OK	OK	-	-	-	-	-	-	-	-	-
SI2C41	OK	OK	-	-	OK	-	OK	-	-	-	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	OK	OK
SL1543	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
SL2325	-	OK	-	-	OK	-	OK	-	-	-	-	OK	OK	-	-	OK	OK	OK	OK	OK	OK	OK	OK
SL3116	-	-	-	-	-	-	-	-	-	-	-	-	OK	-	-	-	-	-	-	-	-	-	-
SL3235	-	OK	-	-	OK	-	OK	-	-	-	-	OK	OK	-	-	OK	OK	OK	OK	OK	OK	OK	OK
SL3323	-	OK	N	N	OK	N	OK	N	N	-	-	OK	OK	-	N	OK	OK	OK	OK	OK	OK	OK	OK
SL3332	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
SP2430	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
SP2442	OK	OK	-	-	OK	-	OK	-	-	-	OK	OK	OK	OK	-	OK	OK	OK	OK	OK	OK	OK	OK
SP2730	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
ST2612	OK	-	-	-	-	-	-	-	-	-	OK	-	OK	OK	-	-	-	-	-	-	-	-	-
ST2613	OK	-	-	-	-	-	-	-	-	-	OK	-	OK	OK	-	-	-	-	-	-	-	-	-
TEE022	OK	N	OK	OK	-	OK	-	-	-	-	N	-	OK	N	OK	-	-	-	-	-	-	OK	-
TLC 3-13 <sup>a</sup>	OK	OK	N	N	OK	N	OK	N	N	-	OK	OK	OK	OK	N	OK	OK	OK	OK	OK	OK	OK	OK
TRT33 <sup>a</sup>	OK	OK	N	N	OK	N	OK	N	N	-	OK	OK	OK	OK	N	OK	OK	OK	OK	OK	OK	OK	OK

a. Whole probe, lumen excluded

Table C-2: Disinfection Agents

	Amuchina	Cavicide	Cavicide1	Cavwipes	Cavwipes1	Cleansept wipes	Mikrobac tissues	Sani Cloth HB wipes	Sani Cloth plus wipes	Sono wipes	Super Sani Cloth wipes	Transeptic	Trophon Companion Cleaning wipes	Umonium medical tissue	Wip'antios Excel
2 MHz CW	-	OK	-	OK	-	-	-	-	-	-	-	-	-	-	-
5 MHz CW	-	OK	-	OK	-	-	-	-	-	-	-	-	-	-	-
AC2541	-	OK	N	OK	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
AL2442	-	OK	N	OK	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
BC441	OK	OK	N	OK	N	OK	OK	OK	OK	OK	N	OK	OK	OK	OK
BE1123	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
BL433	-	-	N	-	N	-	-	-	-	-	N	-	-	-	-
C 1-8	-	OK	N	OK	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
C 2-9	-	N	N	N	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
CA123	OK	OK	N	OK	N	OK	OK	OK	OK	OK	N	OK	OK	OK	OK
CA430E	-	-	N	-	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
CA431	OK	OK	N	OK	N	OK	OK	OK	OK	OK	N	OK	OK	OK	OK
CA541	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
CA631	OK	OK	N	OK	N	OK	OK	OK	OK	OK	N	OK	OK	OK	OK
E 3-12 <sup>a</sup>	-	N	N	N	N	OK	N	OK	N	OK	N	-	OK	OK	OK
EC1123 <sup>a</sup>	-	N	N	N	N	OK	N	OK	N	OK	N	-	OK	OK	OK
EC123 <sup>a</sup>	OK	OK	N	N	N	OK	N	OK	N	OK	N	OK	OK	OK	OK
HF CW	-	OK	-	OK	-	-	-	-	-	-	-	-	-	-	-
IH 6-18	OK	OK	N	OK	N	OK	OK	OK	OK	OK	N	OK	OK	OK	OK
IL 4-13	-	-	N	-	N	-	-	-	-	-	N	-	-	-	-
IOE323	-	-	N	-	N	-	-	-	-	-	N	-	-	-	-
IOT342	-	OK	N	OK	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
L 3-11	-	OK	N	OK	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
L 4-15	-	OK	N	OK	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
L 8-24	-	-	N	-	N	-	-	-	-	-	N	-	-	-	-
LA332E	-	OK	N	OK	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
LA435	-	N	N	N	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
LA523	-	N	N	N	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
LA533	-	OK	N	OK	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
LP323	-	-	N	-	N	-	-	-	-	-	N	-	-	-	-

PROBES AND AGENTS

	Amuchina	Cavicide	Cavicide1	Cavwipes	Cavwipes1	Cleansept wipes	Mikrobac tissues	Sani Cloth HB wipes	Sani Cloth plus wipes	Sono wipes	Super Sani Cloth wipes	Transeptic	Trophon Companion Cleaning wipes	Umonium medical tissue	Wip'antos Excel
LP 4-13	-	-	N	-	N	-	-	-	-	-	N	-	-	-	-
mC 3-11	-	N	N	N	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
P 1-5	-	OK	N	OK	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
P 2-9	-	OK	N	OK	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
P 3-11	-	OK	N	OK	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
P2 3-11	-	OK	N	OK	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
P2 5-13	-	-	N	-	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
PA023E	-	-	N	-	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
PA121E	-	-	N	-	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
PA122E	-	-	N	-	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
PA230E	-	-	N	-	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
PA240	-	OK	N	OK	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
PA250	-	OK	N	OK	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
S2MCW	-	OK	-	OK	-	-	-	-	-	-	-	-	-	-	-
S2MPW	-	OK	-	OK	-	-	-	-	-	-	-	-	-	-	-
S5MCW	-	OK	-	OK	-	-	-	-	-	-	-	-	-	-	-
SB2C41	OK	OK	N	OK	N	OK	OK	OK	OK	OK	N	OK	OK	OK	OK
SB3123	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
SC3123	-	N	N	N	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
SC3421	OK	OK	N	OK	N	OK	OK	OK	OK	OK	N	OK	OK	OK	OK
SE3123 <sup>a</sup>	-	N	N	N	N	OK	N	OK	N	OK	N	-	OK	OK	OK
SE3133 <sup>a</sup>	-	N	N	N	N	OK	N	OK	N	OK	N	-	OK	OK	OK
SHFCW	-	OK	-	OK	-	-	-	-	-	-	-	-	-	-	-
SI2C41	-	OK	N	OK	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
SL1543	-	OK	N	OK	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
SL2325	-	N	N	N	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
SL3116	-	-	N	-	N	-	-	-	-	-	N	-	-	-	-
SL3235	-	N	N	N	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
SL3323	-	N	N	N	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
SL3332	-	OK	N	OK	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
SP2430	-	OK	N	OK	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK

PROBES AND CONSUMABLES



PROBES AND AGENTS

	Amuchina	Cavicide	Cavicide1	Cavwipes	Cavwipes1	Cleansept wipes	Mikrobac tissues	Sani Cloth HB wipes	Sani Cloth plus wipes	Sono wipes	Super Sani Cloth wipes	Transeptic	Trophon Companion Cleaning wipes	Urmonium medical tissue	Wip'antios Excel
SP2442	-	OK	N	OK	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
SP2730	-	OK	N	OK	N	OK	OK	OK	OK	OK	N	-	OK	OK	OK
ST2612	-	-	N	-	N	-	-	-	-	-	N	-	-	-	-
ST2613	-	-	N	-	N	-	-	-	-	-	N	-	-	-	-
TEE022	-	N	N	N	N	-	-	-	-	-	N	-	-	-	-
TLC 3-13 <sup>a</sup>	-	N	N	N	N	OK	OK	OK	N	OK	N	-	OK	OK	OK
TRT33 <sup>a</sup>	-	N	N	N	N	OK	OK	OK	N	OK	N	-	OK	OK	OK

a. Whole probe, lumen excluded

Table C-3: High-Level Disinfection (HLD) Agents

	Antioxyde 1000	Antigermix SI Hypernova Chromos	Cidex	Cidex OPA	Cidex Plus	Compliance	Gigascept AF	Gigascept FF	Metricide 14	Metricide 28	Metricide OPA	Nu-Cidex	Omnicide 14 NS	Omnicide 28 NS	Peracton	Rely+On Perasafe	Revital Ox Resert XL	Sekusept Aktiv	Sterantios	Trophon 2	Trophon EPR	Umonium38	Virkon	Wavicide-01
2 MHz CW	OK	-	-	-	-	OK	N	-	-	-	-	-	-	-	OK	-	OK	-	-	N	N	-	-	-
5 MHz CW	OK	-	-	-	-	OK	N	-	-	-	-	-	-	-	OK	-	OK	-	-	N	N	-	-	-
AC2541	OK	OK	OK	OK	OK	OK	N	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
AL2442	OK	OK	-	OK	-	OK	N	-	-	-	OK	OK	-	-	OK	-	OK	OK	-	-	-	-	-	-
BC441	OK	OK	-	OK	-	OK	N	OK	-	-	OK	OK	-	-	OK	OK	OK	OK	-	OK	OK	OK	-	-
BE1123	OK	N	-	-	-	OK	N	-	-	-	-	-	-	-	OK	-	OK	-	-	OK	OK	-	-	-
BL433	-	N	-	OK	-	-	N	-	-	-	-	-	-	-	-	OK	OK	-	-	OK	OK	-	-	-
C 1-8	OK	OK	OK	OK	OK	OK	N	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
C 2-9	OK	N	-	OK	-	OK	N	-	-	-	OK	OK	-	-	OK	OK	OK	OK	-	OK	OK	OK	-	-
CA123	OK	-	OK	OK	OK	OK	-	OK	OK	-	OK	OK	OK	OK	OK	OK	OK	OK	-	OK	OK	OK	-	OK
CA430E	OK	OK	-	OK	-	OK	N	N	-	-	OK	-	-	-	OK	-	OK	-	-	-	-	OK	OK	-
CA431	OK	-	-	OK	-	OK	N	OK	-	-	OK	OK	-	-	OK	OK	OK	OK	-	OK	OK	OK	-	-
CA541	OK	OK	OK	OK	OK	OK	N	-	OK	-	OK	-	OK	OK	OK	-	OK	-	-	OK	OK	OK	-	OK
CA631	OK	-	OK	OK	OK	OK	N	OK	-	-	OK	OK	OK	OK	OK	OK	OK	OK	-	OK	OK	OK	-	OK
E 3-12 <sup>a</sup>	N	N	-	OK	-	N	N	-	-	-	N	N	-	-	N	OK	OK	N	-	OK	OK	OK	-	-
EC1123 <sup>a</sup>	N	N	-	OK	-	N	N	-	-	-	N	N	-	-	N	OK	OK	N	-	OK	OK	OK	-	-
EC123 <sup>a</sup>	N	N	OK	OK	OK	N	N	N	OK	-	N	N	OK	OK	N	OK	N	N	-	OK	OK	OK	-	OK
HF CW	OK	-	-	-	-	OK	N	-	-	-	OK	-	-	-	OK	-	OK	-	-	N	N	-	-	-
IH 6-18	OK	-	OK	OK	OK	OK	N	OK	OK	OK	OK	OK	-	-	OK	OK	OK	OK	-	-	-	OK	-	-
IL 4-13	-	N	-	-	-	-	N	-	-	-	N	-	-	-	-	-	-	-	N	N	N	-	-	-
IOE323	-	N	-	-	-	-	N	-	-	-	-	-	-	-	-	-	-	-	N	N	N	-	-	-
IOT342	OK	OK	-	OK	-	OK	N	-	-	-	OK	OK	-	-	OK	-	OK	OK	-	-	-	-	-	-
L 3-11	OK	OK	-	OK	-	OK	N	-	-	-	OK	OK	-	-	OK	OK	OK	OK	-	-	-	OK	-	-

PROBES AND CONSUMABLES

C-9

PROBES AND AGENTS

PROBES AND CONSUMABLES

	Anioxyde 1000	Antigermix S1 Hypernova Chronos	Cidex	Cidex OPA	Cidex Plus	Compliance	Gigasept AF	Gigasept FF	Metricide 14	Metricide 28	Metricide OPA	Nu-Cidex	Omnicide 14 NS	Omnicide 28 NS	Peraction	Rely+On Perasafe	Revital Ox Resert XL	Sekusept Aktiv	Sterantos	Trophon 2	Trophon EPR	Uromonium38	Virkon	Wavicide-01
L 4-15	OK	OK	-	OK	-	OK	N	-	-	-	OK	OK	-	-	OK	OK	OK	OK	-	-	-	OK	-	-
L 8-24	-	N	-	OK	-	-	N	-	-	-	-	-	-	-	-	OK	OK	-	-	OK	OK	-	-	-
LA332E	OK	OK	OK	OK	OK	OK	N	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
LA435	-	-	-	-	-	-	N	-	-	-	-	-	-	-	-	-	-	-	-	OK	OK	OK	-	-
LA523	-	OK	-	OK	-	-	N	-	-	-	OK	-	-	-	-	-	-	-	OK	OK	OK	OK	-	-
LA533	OK	OK	OK	OK	OK	OK	N	-	-	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
LP323	-	N	-	-	-	-	N	-	-	-	-	-	-	-	N	-	-	-	N	N	N	-	-	-
LP 4-13	-	N	-	-	-	-	N	-	-	-	-	-	-	-	-	-	-	-	N	N	N	-	-	-
mC 3-11	OK	N	-	OK	-	OK	N	-	-	-	OK	OK	-	-	OK	OK	OK	OK	-	OK	OK	OK	-	-
P 1-5	OK	OK	OK	OK	OK	OK	N	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
P 2-9	OK	OK	-	OK	-	OK	N	-	-	-	OK	OK	-	-	OK	OK	OK	OK	-	-	-	OK	-	-
P 3-11	OK	OK	OK	OK	OK	OK	N	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
P2 3-11	OK	OK	-	OK	-	OK	N	-	-	-	OK	OK	-	-	OK	OK	OK	OK	-	-	-	OK	-	-
P2 5-13	OK	OK	-	OK	-	OK	N	N	-	-	OK	-	-	-	OK	-	OK	-	-	OK	OK	OK	-	-
PA023E	OK	OK	-	OK	-	OK	N	N	-	-	OK	-	-	-	OK	-	OK	-	-	OK	OK	OK	OK	-
PA121E	OK	OK	-	OK	-	OK	N	N	-	-	OK	-	-	-	OK	-	OK	-	-	OK	OK	OK	-	OK
PA122E	OK	OK	-	OK	-	OK	N	N	-	-	OK	-	-	-	OK	-	OK	-	-	OK	OK	OK	OK	-
PA230E	OK	OK	-	OK	-	OK	N	N	-	-	OK	-	-	-	OK	-	OK	-	-	-	-	OK	OK	-
PA240	OK	OK	OK	OK	OK	OK	N	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK
PA250	OK	OK	OK	OK	OK	OK	N	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
S2MCW	OK	-	-	-	-	OK	N	-	-	-	-	-	-	-	OK	-	OK	-	-	N	N	-	-	-
S2MPW	OK	-	-	-	-	OK	N	-	-	-	-	-	-	-	OK	-	OK	-	-	N	N	-	-	-
S5MCW	OK	-	-	-	-	OK	N	-	-	-	-	-	-	-	OK	-	OK	-	-	N	N	-	-	-
SB2C41	OK	OK	-	OK	-	OK	N	OK	-	-	OK	OK	-	-	OK	OK	OK	OK	-	OK	OK	OK	-	-
SB3123	OK	N	-	-	-	OK	N	-	-	-	-	-	-	-	OK	-	OK	-	-	OK	OK	-	-	-

	Anioxyde 1000	Antigermix SI Hypernova Chronos	Cidex	Cidex OPA	Cidex Plus	Compliance	Gigasept AF	Gigasept FF	Metricide 14	Metricide 28	Metricide OPA	Nu-Cidex	Omnicide 14 NS	Omnicide 28 NS	Peraction	Rely+On Perasafe	Revital Ox Resert XL	Sekusept Aktiv	Sterantos	Trophon 2	Trophon EPR	Urmonium38	Virkon	Wavicide-01
SC3123	OK	N	-	OK	-	OK	N	-	-	-	OK	OK	-	-	OK	OK	OK	OK	-	OK	OK	OK	-	-
SC3421	OK	OK	OK	OK	OK	OK	-	OK	-	-	OK	OK	-	-	OK	OK	OK	OK	-	-	-	OK	-	-
SE3123 <sup>a</sup>	N	N	-	OK	-	N	N	-	-	-	N	N	-	-	N	OK	OK	N	-	OK	OK	OK	-	-
SE3133 <sup>a</sup>	N	N	-	OK	-	N	N	-	-	-	N	N	-	-	N	OK	OK	N	-	OK	OK	OK	-	-
SHFCW	OK	-	-	-	-	OK	N	-	-	-	-	-	-	-	OK	-	OK	-	-	N	N	-	-	-
SI2C41	OK	N	-	OK	-	OK	N	-	-	-	OK	OK	-	-	OK	OK	OK	OK	-	-	-	OK	-	-
SL1543	OK	OK	OK	OK	OK	OK	N	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
SL2325	-	-	-	-	-	-	N	-	-	-	-	-	-	-	-	-	-	-	-	OK	OK	OK	-	-
SL3116	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
SL3235	-	-	-	-	-	-	N	-	-	-	-	-	-	-	-	-	-	-	-	OK	OK	OK	-	-
SL3323	-	OK	-	OK	-	-	N	-	-	-	OK	-	-	-	-	-	-	-	OK	OK	OK	OK	-	-
SL3332	OK	OK	OK	OK	OK	OK	N	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
SP2430	OK	OK	OK	OK	OK	OK	N	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	-
SP2442	OK	OK	-	OK	-	OK	N	-	-	-	OK	OK	-	-	OK	OK	OK	OK	-	-	-	OK	-	-
SP2730	OK	OK	OK	OK	OK	OK	N	-	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-	OK
ST2612	OK	N	-	OK	-	-	-	-	-	-	OK	OK	-	-	-	OK	-	OK	OK	N	N	-	-	OK
ST2613	OK	N	-	OK	-	-	-	-	-	-	OK	OK	-	-	-	OK	-	OK	OK	N	N	-	-	OK
TEE022	N	N	OK	OK	OK	-	OK	OK	OK	OK	N	N	OK	OK	-	OK	N	N	OK	N	N	-	OK	OK
TLC 3-13 <sup>a</sup>	N	N	-	OK	-	N	N	-	-	-	N	-	-	-	N	-	OK	-	N	OK	OK	OK	-	-
TRT33 <sup>a</sup>	N	N	-	OK	-	N	N	-	-	-	N	-	-	-	N	-	OK	-	N	OK	OK	OK	-	-

a. Whole probe, lumen included for liquid agents, lumen excluded for Trophon systems. To effectively disinfect lumen use liquid HDL agents following the instruction on chapter 2.

Table C-4: Sterilization Agents

	Steris	Sterrad 100S	Sterrad NX	Sterizone VP4
AC2541	OK	OK	OK	-
AL2442	OK	OK	OK	-
C 1-8	OK	OK	OK	-
E 3-12 <sup>a</sup>	-	-	-	OK
EC1123 <sup>a</sup>	OK	-	OK	OK
EC123 <sup>a</sup>	OK	OK	OK	OK
IH 6-18	OK	-	OK	-
IL 4-13	OK	-	-	-
IOE323	OK	-	-	-
IOT342	OK	OK	OK	-
L 3-11	OK	OK	OK	-
L 4-15	OK	OK	OK	-
LA332E	OK	OK	OK	-
LA523	-	-	-	OK
LA533	OK	OK	OK	-
LP323	OK	-	-	-
LP 4-13	OK	-	-	-
mC 3-11	OK	-	OK	OK
P 1-5	OK	OK	OK	-
P 2-9	OK	OK	OK	-
P 3-11	OK	OK	OK	-
P2 3-11	OK	OK	OK	-
PA240	OK	OK	OK	-
PA250	OK	OK	OK	-
SC3123	OK	-	OK	OK
SE3123 <sup>a</sup>	OK	-	OK	OK
SE3133 <sup>a</sup>	OK	-	OK	OK
SI2C41	OK	OK	OK	-
SL1543	OK	OK	OK	-
SL3116	OK	-	-	-
SL3323	-	-	-	OK
SL3332	OK	OK	OK	-
SP2430	OK	OK	OK	-

	Steris	Sterrad 100S	Sterrad NX	Sterizone VP4
SP2442	OK	OK	OK	-
SP2730	OK	OK	OK	-
TLC 3-13 <sup>a</sup>	-	-	-	OK
TRT33 <sup>a</sup>	-	-	-	OK

a. Whole probe, lumen included

Table C-5: Additional information on agents

Trade name	Type	Chemical type	Manufacturer
Accel RTU wipes	Cleaning	Hydrogen peroxide	Virox Technologies
Amuchina	Disinfection	Sodium Hydroxide	Amuchina / Angelini
Anioxyde 1000	HLD	Hydrogen peroxide	Laboratoires Anios
Antigermix S1 Hypernova Chronos	HLD	UV Radiation	Germitec
Asepti-wipe II	Cleaning	Isopropanol, 2-butoxyethanol, n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chloride	Ecolab
Cavicide	Disinfection	Isopropanol, Ethylene Glycol, NH <sub>4</sub> Cl <sup>a</sup>	Metrex Research Corporation
Cavicide1	Disinfection	Isopropanol, Ethylene Glycol, NH <sub>4</sub> Cl <sup>a</sup> , Ethanol	Metrex Research Corporation
Caviwipes	Disinfection	Isopropanol, Ethylene Glycol, NH <sub>4</sub> Cl <sup>a</sup>	Metrex Research Corporation
Caviwipes1	Disinfection	Isopropanol, Ethylene Glycol, NH <sub>4</sub> Cl <sup>a</sup> , Ethanol	Metrex Research Corporation
Cidex	HLD	Glutaraldehyde	Advanced Sterilization Products (ASP)
Cidex Plus	HLD	Glutaraldehyde	Advanced Sterilization Products (ASP)
Cidex OPA	HLD	Orto-phthalaldehyde	Advanced Sterilization Products (ASP)
Cidezime	Cleaning	Enzimatic detergent	Advanced Sterilization Products (ASP)
CidezimeXTRA	Cleaning	Enzimatic detergent	Advanced Sterilization Products (ASP)
Cleanisept wipes	Disinfection	NH <sub>4</sub> Cl <sup>a</sup>	Dr.Schumacher GmbH
Cleanisept wipes forte	Cleaning	NH <sub>4</sub> Cl <sup>a</sup>	Dr.Schumacher GmbH
Compliance	HLD	Hydrogen peroxide, peracetic acid	Metrex Research Corporation

PROBES AND AGENTS

Trade name	Type	Chemical type	Manufacturer
Enzol	Cleaning	Enzymatic detergent	Advanced Sterilization Products (ASP)
Gigasept AF	HLD	Dimethyldioctyl-ammonium chloride, Phenoxypropanols, Alkylguanidine acetate, Laurylpropylenediamine	Schülke&Mayr GmbH
Gigasept FF	HLD	Reaction product of DMO-THF; ethanol and water	Schülke&Mayr GmbH
Hexanios G+R	Cleaning	NH <sub>4</sub> Cl <sup>a</sup> based	Laboratoires Anios
Klenzyme	Cleaning	Enzymatic detergent	Steris Corporation
Metricide 14	HLD	Glutaraldehyde	Metrex Research Corporation
Metricide 28	HLD	Glutaraldehyde	Metrex Research Corporation
Metricide OPA	HLD	Orto-phtalaldehyde	Metrex Research Corporation
Metricide Plus 30	HLD	Glutaraldehyde	Metrex Research Corporation
Metrizyme	Cleaning	Enzymatic detergent	Metrex Research Corporation
Microzid AF wipes	Cleaning	Ethanol, Propan-1-ol	Schülke&Mayr GmbH
Microzid PAA wipes	Cleaning	Hydrogen peroxide, peracetic acid	Schülke&Mayr GmbH
Microzid Sensitive wipes	Cleaning	NH <sub>4</sub> Cl <sup>a</sup>	Schülke&Mayr GmbH
Mild soap	Cleaning		Various
Nu-Cidex	HLD	Orto-phtalaldehyde	Advanced Sterilization Products (ASP)
Omnicide 14 NS	HLD	Glutaraldehyde	Metrex Research Corporation
Omnicide 28 NS	HLD	Glutaraldehyde	Metrex Research Corporation
Peraction	HLD	Hydrogen peroxide	Biosteril
Percept RTU cleaner	Cleaning	Hydrogen peroxide	Virox Technologies
Prolystica	Cleaning	Enzymatic detergent	Steris Corporation
Protex wipes	Disinfection	NH <sub>4</sub> Cl <sup>a</sup>	Parker Laboratories Inc.
Rely+On Perasafe	HLD	Sodium Percarbonate, Citric acid, Sodium Carbonate	Chemours
Revital Ox Resert XL	HLD	Hydrogen peroxide	Steris Corporation
Sani-Cloth Active wipes	Cleaning	NH <sub>4</sub> Cl <sup>a</sup>	Nice-Pak / Professional Disposables International (PDI)
Sani Cloth HB wipes	Disinfection	NH <sub>4</sub> Cl <sup>a</sup>	Nice-Pak / Professional Disposables International (PDI)
Sani Cloth plus wipes	Disinfection	NH <sub>4</sub> Cl <sup>a</sup> , Isopropanol	Nice-Pak / Professional Disposables International (PDI)
SaniZide plus (spray, wipes)	Cleaning	NH <sub>4</sub> Cl <sup>a</sup>	Safetec

PROBES AND AGENTS

Trade name	Type	Chemical type	Manufacturer
Sono wipes	Disinfection	NH <sub>4</sub> Cl <sup>a</sup>	Advanced Ultrasound Solutions
Sekusept Aktiv	HLD	Sodium Percarbonate, Citric acid, Sodium Carbonate	Ecolab
Steranios	HLD	Glutaraldehyde	Laboratoires Anios
Steris	Sterilization	Hydrogen peroxide	Steris Corporation
Sterrad 100S	Sterilization	Hydrogen peroxide	Advanced Sterilization Products (ASP)
Sterrad NX	Sterilization	Hydrogen peroxide	Advanced Sterilization Products (ASP)
Sterizone VP4	Sterilization	Hydrogen peroxide, Ozone	TSO <sub>3</sub>
Super Sani Cloth wipes	Disinfection	NH <sub>4</sub> Cl <sup>a</sup> , Isopropyl Alcohol	Nice-Pak / Professional Disposables International (PDI)
Transeptic	Disinfection	Isopropyl alcohol	Parker Laboratories Inc.
Trionic D (spray, wipes)	Cleaning	NH <sub>4</sub> Cl <sup>a</sup>	Ebiox
Tristel sporicidal wipes	Cleaning	Chlorine dioxide	Tristel Solution Ltd.
Trophon 2	HLD	Hydrogen peroxide	Nanosonic Ltd.
Trophon Companion Cleaning wipes	Disinfection	NH <sub>4</sub> Cl <sup>a</sup>	Nanosonic Ltd.
Trophon EPR	HLD	Hydrogen peroxide	Nanosonic Ltd.
Umonium38	HLD	NH <sub>4</sub> Cl <sup>a</sup>	Laboratoire Huckert's International
Umonium medical tissue	Disinfection	NH <sub>4</sub> Cl <sup>a</sup>	Laboratoire Huckert's International
Virusolve + wipes	Cleaning	2-Aminoethanol, NH <sub>4</sub> Cl <sup>a</sup>	Amity International
Virkon	HLD	Pentapotassium bis(peroxymonosulphate) bis(sulphate)	Day-Impex Ltd
Wavicide-01	HLD	Glutaraldehyde	Wave Energy Systems
Wip'anios Excel	Disinfection	NH <sub>4</sub> Cl <sup>a</sup>	Laboratoires Anios

a. Ammonium Chloride