

Connect your MyLab to your PC, easily!

MyLab™ Desk



- > Archive, review, post-process, print exams from a PC
- > All the comfort, functions and convenience of the PC
- > MyLab-based user interface

MyLab Desk

Image management made easy


MyLab Desk is Esaote's answer to its users' need for a simple, straightforward way to archive, review, post-process, report and print their MyLab examinations on a PC in the comfort of their (home) office or while travelling.

MyLab Desk is fully integrated into Esaote's MyLab platform. The program can be easily downloaded from the ultrasound system to any Windows® XP (Home/Professional)-based PC or laptop via the USB personal memory drive or CD/DVD. MyLab Desk's user interface and workflow is very similar to the MyLab systems', so the program is easy-to-learn and use for MyLab users. MyLab Desk provides a way to increase workflow and productivity in private offices, as well as in clinic and hospital departments.

MyLab Desk makes it possible to:

- > Archive, review and post-process examinations performed with the MyLab ultrasound systems.
- > Import native Esaote file formats (raw data UAF and EAF) via USB, CD/DVD and network.
- > Perform generic and application-specific measurements.
- > Review, modify and print the examinations (reports and images).
- > Export data using the standard features of the PC, i.e. burn on a CD/DVD, email, etc.

MyLab Desk offers the modern ultrasound professional a flexible way to connect their MyLab to their PC, easily!



833 01719 000 (MA Rev. B)