



IHE Integration Statement		Date: Oct. 26, 2007
Vendor	Product Name	Version
ESAOTE S.p.A.	MyLab 25/30/40/50/50 XV * * all with DICOM® option	5.01 / 6.01 / 6.02 / 6.03 / 6.10 / 6.11 / 6.12 / 7.01 / 7.02 / 7.10 / 7.20
	MyLab 60/70/70 XVG * * all with DICOM® option	1.90 / 2.00 / 2.01 / 2.02 / 2.03 / 2.04 / 2.05 / 3.00 / 3.01 / 3.10 / 3.11
This product implements all of the transactions specified in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	Acquisition Modality	Broad Worklist Query
Patient Information Reconciliation	Acquisition Modality	none
Portable Data for Imaging	Portable Media Creator	none
Echocardiography Workflow	Acquisition Modality	none
Internet address for vendor's IHE information: www.esaote.com/products/ultrasound/statements/IHE2.htm Vendor's e-mail address for information about IHE: ihe@esaote.com		
Links to Standards Conformance Statements for the Implementation		
DICOM®	www.esaote.com/products/ultrasound/statements/confstat.htm	
Links to general information on IHE		
In North America: www.ihe.net		In Europe: www.ihe-europe.org

ESAOTE S.p.A.
via Siffredi 58
16153 Genova (Italy)

© Copyright ESAOTE S.p.A. 2007. All rights reserved.

[®] DICOM is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information.