



| <b><i>IHE Integration Statement – Interoperability Statement</i></b> |  |                       |
|--|--|-----------------------|
| <b><i>Manufacturer</i></b>   | <b><i>Product Name</i></b>                                 | <b><i>Version</i></b> |
| ESAOTE S.p.A.  | SUITESTENSA  | 1.6                   |
| <b><i>Integration Profiles supported</i></b>                         |  |                       |
|  | <b><i>Supported Actors</i></b>                             |                       |
| Scheduled Workflow   | DSS/Order Filler<br>Image Display<br>Image Manager/Archive |                       |
| Patient Information Reconciliation                                   | DSS/Order Filler<br>Image Display<br>Image Manager/Archive |                       |
| Cardiac Catheterization Workflow                                     | Image Display<br>Image Manager/Archive<br>Evidence Creator |                       |
| Echo Workflow  | Image Display<br>Evidence Creator<br>Image Manager/Archive |                       |
| Stress Testing Workflow  | Image Display<br>Evidence Creator<br>Image Manager/Archive |                       |
| Consistent Presentation of Images                                    | Image Display<br>Print Composer<br>Image Manager/Archive   |                       |
| Presentation of Grouped Procedures                                   | Image Manager/Archive                                      |                       |

|                                    |  |
|------------------------------------|--|
| Evidence Documents                 | Image Manager/Archive  |
| Evidence Documents -<br>Cardiology | Image Manager/Archive<br>Image Display                                       |
| Key Image Note                     | Image Manager/Archive<br>Image Display                                       |
| Mammography Image                  | Image Display<br>Image Manager/Archive                                       |
| Portable Data for Imaging          | Report Reader<br>Portable Media Creator<br>Image Display                     |
| XDS.b-i                            | Imaging Document Consumer<br>Imaging Document Source                         |
| ARI                                | Image Display<br>Image Manager/Archive<br>Report Reader<br>Report Repository |
| SINR                               | Report Reader<br>Report Repository   |
| ATNA                               | Secure Application   |
| Consistent Time                    | Time Client  |
| <b><i>Date of Declaration</i></b>  | Dicember 2012  |

**Esote S.p.A. – EBITAET Business Unit**  
**Product Planning and Integration- Director**  
Via Siffredi 58 I-16153 Genova