



Suggested Printers and DVRs

Model	Type	Dimension W,H,D, [mm]	Weight [Kg]	Interfaces
Sony UP-X898MD	Thermal B/W	154 x 88 x 240	2.5	Composite, USB
HP Color LaserJet Pro M454dw	Laser color	412 x 469 x 295	18.9	LAN, Wi-Fi, USB
Xerox C230 Colour Laser	Laser color	368 x 394 x 252	16.1	LAN, Wi-Fi, USB

Please note that the Ultrasound Peripherals are system/configuration dependent.

Please check with your Service representative the compatibility between your MyLab™ system and the desired printer or DVR, and ask for the installation instructions and certified drivers to use.

Installing with wrong drivers could compromise the correct system functionality.

Technology and features are device/configuration-dependent. Specifications subject to change without notice. Information might refer to products or modalities not yet approved in all countries. For further details, please contact your Esaote sales representative.