



evo^{'16}lution

Dedicated MRI




Head imaging, the logical
evolution of Esaote MRI


**Evolution,
Let's Think A-Head**

esaote


The coil design offers a comfortable open view to reduce claustrophobic response. The sliding mechanism makes for easy and fast patient positioning. The flexible Esaote user interface allows for full customization of the standard head protocols. All coils are supplied with paddings for stable and comfortable positioning.




SE T1 Sagittal




FSE T2 Transverse




SE T1 Coronal



FAST FLAIR Sagittal



GE T2 FC Transverse



3D SST1 0,4 mm Transverse

"We have drawn on our thirty years of experience in MSK MRI to enable to our customers to take a step further in the use of our MRI systems"

EVO'16 introduces the head imaging option for both the S-scan and G-scan Brio MRI systems, completing the functionality of both scanners.

EVO'16 makes it possible to extend the clinical application, consequently enlarging the patient database and improving system economics.

The Head Imaging Option is available for new systems and as an upgrade for the installed base (contact local distributor for details).



The G-scan Brio allows for weight-bearing imaging of the head, opening a new field of application in MRI. For weight-bearing imaging, the ergonomic design of the G-scan Brio ensures a comfortable sitting position for the patient.





evo¹⁶lution

Head imaging for G-scan Brio and S-scan



ESAOTE S.p.A.- sole-shareholder company
Via Enrico Melen, 77 16152 Genova, ITALY, Tel. +39 010 6547 1, Fax +39 010 6547 275, info@esaote.com

Technology and features are system/configuration dependent. Specifications subject to change without notice.
Information might refer to products or modalities not yet approved in all countries. Product images are for illustrative purposes only.
For further details, please contact your Esaote sales representative.



Please visit us online
for more information

www.esaote.com