



Never stop striving for  
**e**asy Scanning  
**e**asy Quantification  
in **Critical Care**



N E V E R   S T O P   S E E I N G   T H E   U N S E E N .




# Esaote solutions: make your **critical care examination easier**



## **AutoEF**

Auto LV border tracking in less than 5 seconds to get Simpson Bi-plane EF measurement.



## **@Doppler**

One click to get automatic box positioning, box and line steering, sample positioning and angle correction to facilitate the Doppler acquisition.




## **@Scan**

Continuous automatic gain adjustment to optimize image quality in real time through a smart algorithm, reducing manipulations.



## **microV**

Maximum sensitivity to highlight even the smallest flows and the tissue perfusion of the vascularized part.



## **XStrain™**

Immediate layout of GLS with zero-click technology for clear representation of the segments' contractility and bull's eye score.



## **MyLibrary**

On-board reference guide with anatomical, ecographic, how to scan, and textual description.



Esaote S.p.A. - sole-shareholder company - Via Enrico Meloni 77, 16152 Genova, ITALY, Tel. +39 010 6547 1, Fax +39 010 6547 275, info@esaote.com

Technology and features are system/configuration dependent. Specifications subject to change without notice. Information might refer to products or modalities not yet approved in all countries. Product images are for illustrative purposes only. For further details, please contact your Esaote sales representative.

N E V E R   S T O P   S E E I N G   T H E   U N S E E N .

**esaote**