

Focus on:

EuroEcho2013
Imaging

11-14 December

Istanbul - Turkey

1978-2013

35 Years in Cardiovascular Ultrasound



2003 MYLAB WORLD



TODAY MYLABALPHA



1995 AU3 PARTNER



1999 CARIS/MEGAS



1990 SIM 7000 CFM

1988 SIM 5000 PLUS



1985 SIM 4000



1984 SIM 3000



1980 ECHOPROCESSOR



1982 NEW SECTOR IMAGER





1978 ECHOVISOGRAPH





- Heart-vessel combined examination
- Early CVD detection
- Treatment monitoring

Improved quality of life & more affordable care




PREVENTION Suite



XStrainTM 4D

Get the Whole Picture



XStrainTM is a non-invasive tool to better investigate myocardial function and explore and quantify aspects of the heart's physiology which were not possible to detect and quantify with previous ultrasound technologies.



Get the New DVD!

Click on
esaote.com/education
or e-mail to
info@esaote.com

Focus on:

EuroEcho2013 Imaging

Echocardiography





The top quality ultrasound cardiovascular solution

- Specific cine tools for ultrasounds: different screen tiles, zoom, trim, customizable templates, side-by-side comparison of images or loops
- Dedicated DICOM SR reporting module based on SIEC/ASE classification (standard dataset for normal/pathological exams)
- Bookmark saving of important findings
- Specific datasheets for all pathologies
- Complete 2D quantitative analysis (B-Mode, M-Mode; doppler; calibration) and qualitative analysis (structured morphological descriptions and free text descriptions)
- Advanced functionalities including complete TTE, TEE
- Analysis of technical raw data, in real-time, based on customizable criteria
- Vendor-independent advanced 2D,3D & 4D processing and reconstruction ultrasound analysis packages
- Automatic PDF report generation, with patient data, key images, clinical and quantitative data, reconstructions



CVIS PACS



iQProbe transducer technology

iQProbes

iQProbe transducer technology is based on the following advanced developments:

- **New active PZT material and Resin**
- **5 adaptive matching layers**
- **New filling material**
- **New acoustic lens**

Focus on:

EuroEcho2013 Imaging

11-14 December
Istanbul - Turkey

Esaote, your partner in Cardiology

MyLabTwice



MyLabSat

MyLabSeven



MyLabAlpha



MyLab40



MyLabOne



MyCardioPad



Come and **play** with us
& **win** a free course!

International advanced course on
Echocardiography Deformation Indexes

Scientific organiser

**Dr. Claudio
Bussadori**

Department of Pediatric
Cardiology and Adults
with Congenital Heart
Disease, IRCCS San
Donato Hospital Milano,
Milan, Italy



Committed to Education

Scientific secretariat

Serono Symposia International Foundation Salita di San Nicola da Tolentino, 1/b
00187 Rome, Italy - Phone: +39-06-420413 302 Fax: +39-06-420413 677
E-mail: info@seronosymposia.org

Serono Symposia International Foundation is a Swiss Foundation with headquarters in 14,
rue du Rhône, 1204 Geneva, Switzerland



Esaote S.p.A.

Via di Caciolle, 15 50127 Florence, Italy, Tel. +39 055 4229 1, Fax +39 055 4229 208, international.sales@esaote.com
Via A. Siffredi, 58 16153 Genoa, Italy, Tel. +39 010 6547 1, Fax +39 010 6547 275, info@esaote.com

FRANCE
Esaote France S.A.R.L.
ZA du Bel Air - 10, rue de Témara
78105 Saint-Germain-en-Laye
Tel. +33 1 8204 8900
Fax +33 1 3061 7210
info@esaote.fr

BRASIL
Brazilian Direct Office
Rua Tomas Carvalhal, 711
04006-001 São Paulo SP
Tel. +55 11 2589 0533
Fax +55 11 2589 0527
leonardo.pili@esaote.com.br

GERMANY
Esaote Biomedica Deutschland GmbH
Max-Planck-Straße 27a
50858 Köln
Tel. +49 2234 688 5600
Fax +49 2234 967 9628
info@esaote.de

ARGENTINA
Esaote Latinoamérica S.A.
San Martín 551, Cuerpo 'C', Piso 8,
(C1004AAK) Buenos Aires
Tel. +54 11 4326 1832
Fax: +54 11 4328 1245
info@esaote.com.ar

SPAIN
Esaote España S.A.
C/. Pont Reixat, 5
08960 Sant Just Desvern, Barcelona
Tel. +34 93 473 2090
Fax +34 93 473 2042
info@esaote.es

INDIA
Esaote Asia Pacific Diagnostic Private Limited
DLF IT Park, A - 44 & 45, Tower- C,
Ground Floor, Sector-62, Noida,
Uttar Pradesh, India
Pin Code: 201 301, Tel. +91 120 4732444
Fax +91 120 4750148
info@esaote.in

THE NETHERLANDS AND BELGIUM
Esaote Benelux B.V.
Philipsweg 1
6227 AJ Maastricht
Tel. +31 43 3824650
Fax +31 43 3824651
benelux@esaote.nl

HONG KONG AND FAR EAST
Esaote China Ltd
18/F, 135 Bonham Strand Trade Centre,
135 Bonham Strand, Sheung Wan, Hong Kong
Tel. +852 2545 8386
Fax +852 2543 3068
esaote@esaotechina.co

UK
Esaote UK
14, Cambridge Science Park
Milton Road, Cambridge, CB4 0FQ
Tel. +44 1223 424499
Fax +44 709 288 0231
infoUK@esaote.com

CHINA
Esaote Shenzhen Medical Equipment
Room 2608, Tower B
Beijing Global Trade Center
36 North Third Ring Road East,
Dongcheng District, 100013, Beijing
Tel. +86 010 58257766
Fax +86 010 52257760

NORTH AMERICA
Esaote North America
8000 Castleway Drive,
Indianapolis, IN 46250
Tel. +1 317 813 6000
Fax +1 317 813 6600
inquire@esaoteusa.com

RUSSIAN FEDERATION AND CIS
Esaote S.p.A.
18 Leningradsky prospekt
Off. 5 and 6, Moscow 125040
Tel. +7 495 232 0205
Fax +7 495 232 1833
esaotemoscw@yandex.ru

www.esaote.com